

For Research Use Only

SCD Recombinant antibody, PBS Only (Detector)

Catalog Number: 81468-3-PBS



Basic Information

Catalog Number: 81468-3-PBS	GenBank Accession Number: BC005807	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6319	CloneNo.: 241487E11
Source: Rabbit	UNIPROT ID: O00767	
Isotype: IgG	Full Name: stearoyl-CoA desaturase (delta-9-desaturase)	
Immunogen Catalog Number: AG29156	Calculated MW: 355 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01272-1, SCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81468-2-PBS. Detection antibody: 81468-3-PBS. Standard: Ag29156. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01272-2, SCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81468-4-PBS. Detection antibody: 81468-3-PBS. Standard: Ag29156. Range: 0.313-40 ng/mL