

For Research Use Only

DNMT3A Recombinant antibody, PBS Only (Detector)

Catalog Number: 81474-4-PBS



Basic Information

Catalog Number: 81474-4-PBS	GenBank Accession Number: BC043617	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1788	CloneNo.: 242257F5
Source: Rabbit	UNIPROT ID: Q9Y6K1	
Isotype: IgG	Full Name: DNA (cytosine-5-)-methyltransferase	
Immunogen Catalog Number: AG15111	3 alpha	
	Calculated MW: 912 aa, 102 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01598-2, DNMT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81474-2-PBS. Detection antibody: 81474-4-PBS. Standard: Ag15111. Range: 2.5-80 ng/mL.