

For Research Use Only

PPARG Recombinant antibody, PBS Only (Detector)



Catalog Number: 81490-4-PBS

Basic Information

Catalog Number: 81490-4-PBS	GenBank Accession Number: BC006811	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5468	CloneNo.: 230374E7
Source: Rabbit	UNIPROT ID: P37231	
Isotype: IgG	Full Name: peroxisome proliferator-activated receptor gamma	
Immunogen Catalog Number: AG10005	Calculated MW: 58 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00256-1, PPAR γ Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81490-3-PBS. Detection antibody: 81490-4-PBS. Standard: Ag10005. Range: 1.25-20 ng/mL