

For Research Use Only

PPAR Gamma Recombinant antibody, PBS Only (Capture)

Catalog Number: 81490-7-PBS



Basic Information

Catalog Number: 81490-7-PBS	GenBank Accession Number: BC006811	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5468	CloneNo.: 230374C4
Source: Rabbit	UNIPROT ID: P37231	
Isotype: IgG	Full Name: peroxisome proliferator-activated receptor gamma	
Immunogen Catalog Number: AG10005	Calculated MW: 58 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00256-3, Human PPAR Gamma Recombinant Matched Antibody Pair - PBS only. 81490-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10005. 81490-6-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



Serum of sixteen individual healthy human donors was measured. The PPAR Gamma concentration of detected samples was determined to be 7.05 ng/mL with a range of 0.29 - 35.09 ng/mL.