

For Research Use Only

# SNAI1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 81584-2-PBS

## Basic Information

Catalog Number:

81584-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24248

GenBank Accession Number:

BC012910

GeneID (NCBI):

6615

UNIPROT ID:

O95863

Full Name:

snail homolog 1 (Drosophila)

Calculated MW:

264 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

240164A5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00343-1, SNAI1 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 81584-2-PBS. Detection antibody: 81584-3-PBS. Standard: Ag24248. Range: 0.313-20 ng/mL.