

For Research Use Only

GAPDH Recombinant antibody, PBS Only (Capture)



Catalog Number: 81640-3-PBS

Basic Information

Catalog Number: 81640-3-PBS	GenBank Accession Number: BC004109	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2597	CloneNo.: 2L15
Source: Rabbit	UNIPROT ID: P04406	
Isotype: IgG	Full Name: glyceraldehyde-3-phosphate dehydrogenase	
Immunogen Catalog Number: AG0766	Calculated MW: 36 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00549-1, GAPDH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81640-3-PBS. Detection antibody: 81640-5-PBS. Standard: ag0766. Range: 0.098-100 ng/mL.