

For Research Use Only

CHMP4B Recombinant antibody, PBS Only (Detector)



Catalog Number: 81691-4-PBS

Basic Information

| | | |
|--|---|---|
| Catalog Number: 81691-4-PBS | GenBank Accession Number: BC033859 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 128866 | CloneNo.: 240857E4 |
| Source: Rabbit | UNIPROT ID: Q9H444 | |
| Isotype: IgG | Full Name: chromatin modifying protein 4B | |
| Immunogen Catalog Number: AG4544 | Calculated MW: 224 aa, 25 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00733-1, CHMP4B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81691-3-PBS. Detection antibody: 81691-4-PBS. Standard: Ag4544. Range: 1.25-80 ng/mL.