

For Research Use Only

BCL6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 81751-1-PBS

Basic Information

Catalog Number: 81751-1-PBS	GenBank Accession Number: BC150184	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 604	CloneNo.: 5N11
Source: Rabbit	UNIPROT ID: P41182	
Isotype: IgG	Full Name: B-cell CLL/lymphoma 6	
Immunogen Catalog Number: AG15519	Calculated MW: 706 aa, 79 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80008-1, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-1-PBS. Detection antibody: 81751-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL