

For Research Use Only

4EBP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 81812-3-PBS

Basic Information

Catalog Number: 81812-3-PBS	GenBank Accession Number: BC004459	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1978	CloneNo.: 240535G7
Source: Rabbit	UNIPROT ID: Q13541	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 4E binding protein 1	
	Calculated MW: 118 aa, 12 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00559-1, Phospho-4EBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81812-2-PBS. Detection antibody: 81812-3-PBS. Standard: SY01767P. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP00559-2, Phospho-4EBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81812-4-PBS. Detection antibody: 81812-3-PBS. Standard: SY01767P. Range: 0.78-50 ng/mL.