

For Research Use Only

RAB1A Recombinant antibody, PBS Only



Catalog Number: 81837-1-PBS

Basic Information

Catalog Number:

81837-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2264

GenBank Accession Number:

BC000905

GeneID (NCBI):

5861

UNIPROT ID:

P62820

Full Name:

RAB1A, member RAS oncogene family

Calculated MW:

205 aa, 23 kDa

Observed MW:

23-25 kDa

Purification Method:

Protein A purification

CloneNo.:

2G24

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81837-1-RR (RAB1A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81837-1-PBS in a different storage buffer formulation.