

For Research Use Only

ANGPT1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 81990-2-PBS

Basic Information

Catalog Number: 81990-2-PBS	GenBank Accession Number: BC152411	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 284	CloneNo.: 240764A11
Source: Rabbit	UNIPROT ID: Q15389	
Isotype: IgG	Full Name: angiopoietin 1	
Immunogen Catalog Number: AG18829	Calculated MW: 498 aa, 58 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00754-1, Angiopoietin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81990-2-PBS. Detection antibody: 81990-3-PBS. Standard: Ag18829. Range: 1.56-100 ng/mL.