

For Research Use Only

PINK1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 81991-3-PBS

Basic Information

Catalog Number: 81991-3-PBS	GenBank Accession Number: BC028215	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 65018	CloneNo.: 240449F7
Source: Rabbit	UNIPROT ID: Q9BXM7	
Isotype: IgG	Full Name: PTEN induced putative kinase 1	
Immunogen Catalog Number: AG19825	Calculated MW: 581 aa, 63 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00404-2, PINK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81991-5-PBS. Detection antibody: 81991-3-PBS. Standard: Ag19825. Range: 0.78-100 ng/mL