

For Research Use Only

CHAT Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 81994-4-PBS

Basic Information

Catalog Number:

81994-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_020549

GeneID (NCBI):

1103

UNIPROT ID:

P28329

Full Name:

choline acetyltransferase

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

240604F5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00576-1, CHAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81994-4-PBS. Detection antibody: 81994-3-PBS. Standard: SY00995. Range: 19.5-1250 pg/mL



Cytometric bead array standard curve of MP00576-2, CHAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81994-2-PBS. Detection antibody: 81994-4-PBS. Standard: SY00995. Range: 19.5-2500 pg/mL