

For Research Use Only

G3BP2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82080-3-PBS

Basic Information

Catalog Number: 82080-3-PBS	GenBank Accession Number: BC011731	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9908	CloneNo.: 230275B3
Source: Rabbit	UNIPROT ID: Q9UN86	
Isotype: IgG	Full Name: GTPase activating protein (SH3 domain) binding protein 2	
Immunogen Catalog Number: AG9355	Calculated MW: 482aa,54 kDa; 449aa,51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00090-1, G3BP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82080-4-PBS. Detection antibody: 82080-3-PBS. Standard: Ag9355. Range: 0.625-20 ng/mL