

For Research Use Only

# CD44 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82188-3-PBS

## Basic Information

<b>Catalog Number:</b> 82188-3-PBS	<b>GenBank Accession Number:</b> BC004372	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 960	<b>CloneNo.:</b> 230326G3
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000026508	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P16070	
<b>Immunogen Catalog Number:</b> AG8143	<b>Full Name:</b> CD44 molecule (Indian blood group)	
	<b>Calculated MW:</b> 742 aa, 82 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00033-1, CD44  
Recombinant Matched Antibody Pair, PBS Only.  
Capture antibody: 82188-3-PBS. Detection  
antibody: 82188-4-PBS. Standard: Ag8143. Range:  
78-5000 pg/mL.