

For Research Use Only

CD44 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82188-4-PBS

Basic Information

Catalog Number: 82188-4-PBS	GenBank Accession Number: BC004372	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 960	CloneNo.: 230326G9
Source: Rabbit	ENSEMBL Gene ID: ENSG00000026508	
Isotype: IgG	UNIPROT ID: P16070	
Immunogen Catalog Number: AG8143	Full Name: CD44 molecule (Indian blood group)	
	Calculated MW: 742 aa, 82 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00033-1, CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82188-3-PBS. Detection antibody: 82188-4-PBS. Standard: Ag8143. Range: 78-5000 pg/mL



Sandwich ELISA standard curve of MP00033-2, CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82188-2-PBS. Detection antibody: 82188-4-PBS. Standard: Ag8143. Range: 78-5000 pg/mL