

For Research Use Only

Aquaporin 4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82296-4-PBS



Basic Information

Catalog Number: 82296-4-PBS	GenBank Accession Number: BC022286	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 361	CloneNo.: 240856E10
Source: Rabbit	UNIPROT ID: P55087	
Isotype: IgG	Full Name: Aquaporin 4	
Immunogen Catalog Number: AG9561	Calculated MW: 323 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00781-3, AQP4/Aquaporin 4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82296-4-PBS. Detection antibody: 82296-3-PBS. Standard: Ag9561. Range: 0.313-40 ng/mL