

For Research Use Only

Cytokeratin 5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82548-1-PBS

Basic Information

Catalog Number: 82548-1-PBS	GenBank Accession Number: BC024292	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3852	CloneNo.: 2I12
Source: Rabbit	UNIPROT ID: P13647	
Isotype: IgG	Full Name: keratin 5	
Immunogen Catalog Number: AG24184	Calculated MW: 590 aa, 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00650-1, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-18-PBS. Detection antibody: 82548-1-PBS. Standard: Ag24184. Range: 0.625-80 ng/mL