

For Research Use Only

# KRT5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82548-8-PBS

## Basic Information

Catalog Number: 82548-8-PBS	GenBank Accession Number: BC024292	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3852	CloneNo.: 4C2
Source: Rabbit	UNIPROT ID: P13647	
Isotype: IgG	Full Name: keratin 5	
Immunogen Catalog Number: AG24184	Calculated MW: 590 aa, 62 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00650-2, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-21-PBS. Detection antibody: 82548-8-PBS. Standard: Ag24184. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP00650-3, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-4-PBS. Detection antibody: 82548-8-PBS. Standard: Ag24184. Range: 0.625-80 ng/mL