

For Research Use Only

Rat IL2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82626-2-PBS

Basic Information

Catalog Number:

82626-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_053836

GeneID (NCBI):

116562

UNIPROT ID:

P17108

Full Name:

interleukin 2

Calculated MW:

18kd

Purification Method:

Protein A purification

CloneNo.:

240264B2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Rat

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00535-1, RAT IL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82626-2-PBS. Detection antibody: 82626-3-PBS. Standard: Eg0260. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP00535-2, RAT IL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82626-2-PBS. Detection antibody: 82626-4-PBS. Standard: Eg0260. Range: 0.156-20 ng/mL.