

For Research Use Only

VDAC3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82666-14-PBS

Basic Information

Catalog Number:

82666-14-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005662

GeneID (NCBI):

7419

UNIPROT ID:

Q9Y277

Full Name:

voltage-dependent anion channel 3

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

230451C7

Applications

Tested Applications:

WB, Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82666-14-RR (VDAC3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82666-14-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00287-1, VDAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82666-14-PBS. Detection antibody: 82666-15-PBS. Standard: SY01307. Range: 0.156-10 ng/mL.