

For Research Use Only

# VDAC3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82666-15-PBS

## Basic Information

Catalog Number:  
82666-15-PBS

Size:  
1mg/ml

Source:  
Rabbit

Isotype:  
IgG

GenBank Accession Number:  
NM\_005662

GeneID (NCBI):  
7419

UNIPROT ID:  
Q9Y277

Full Name:  
voltage-dependent anion channel 3

Calculated MW:  
31 kDa

Purification Method:  
Protein A purification

CloneNo.:  
230451H4

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00287-1, VDAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82666-14-PBS. Detection antibody: 82666-15-PBS. Standard: SY01307. Range: 0.156-10 ng/mL



Cytometric bead array standard curve of MP00287-2, VDAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82666-15-PBS. Detection antibody: 82666-18-PBS. Standard: SY01307. Range: 0.156-10 ng/mL