

For Research Use Only

# VDAC3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82666-16-PBS

## Basic Information

Catalog Number:

82666-16-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005662

GeneID (NCBI):

7419

UNIPROT ID:

Q9Y277

Full Name:

voltage-dependent anion channel 3

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

230451D4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00287-3, VDAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82666-16-PBS. Detection antibody: 82666-17-PBS. Standard: SY01307. Range: 0.156-10 ng/mL