

For Research Use Only

# MFN2 Recombinant antibody, PBS Only



Catalog Number: 82673-2-PBS

## Basic Information

<b>Catalog Number:</b> 82673-2-PBS	<b>GenBank Accession Number:</b> BC017061	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 9927	<b>CloneNo.:</b> 2A5
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000116688	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> O95140	
<b>Immunogen Catalog Number:</b> AG2845	<b>Full Name:</b> mitofusin 2	
	<b>Calculated MW:</b> 757 aa, 86 kDa	
	<b>Observed MW:</b> 86 kDa	

## Applications

**Tested Applications:**  
WB, IF, IHC, Indirect ELISA

**Species Specificity:**  
Human, rat

## Background Information

Mitofusin-2 (MFN2), a mitochondrial membrane protein, is a dynamin-like GTPase embedded in the outer mitochondrial membrane (OMM) which plays a central role in regulating mitochondrial fusion and cell metabolism (PMID: 35399520; 34026850). There are two human mitofusin isoforms MFN1 and MFN2, MFN1 and MFN2 mediate outer membrane fusion, OPA1 is involved in inner membrane fusion, and DRP1 is responsible for mitochondrial fission (PMID: 17084350; 35399520). It has been shown that MFN2 participates in multiple cellular biological processes under physical and pathological conditions. Apart from mitochondria fusion, MFN2 is also involved in mitochondrial-ER contacts (MERCs), mitophagy, cellular apoptosis, and proliferation (PMID: 34026850).

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 82673-2-RR (MFN2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using MFN2 antibody (82673-2-RR, Clone: 2A5) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 82673-2-RR (MFN2 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 82673-2-RR (MFN2 antibody) at dilution of 1:2000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.