

For Research Use Only

# DNMT1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82678-3-PBS

## Basic Information

<b>Catalog Number:</b> 82678-3-PBS	<b>GenBank Accession Number:</b> BC126227	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 1786	<b>CloneNo.:</b> 240302H8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P26358	
<b>Isotype:</b> IgG	<b>Full Name:</b> DNA (cytosine-5-)-methyltransferase	
<b>Immunogen Catalog Number:</b> AG29479	<b>Calculated MW:</b> 1632 aa, 185 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00377-1, DNMT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82678-3-PBS. Detection antibody: 82678-2-PBS. Standard: Ag29479. Range: 1.56-100 ng/mL