

For Research Use Only

# HRAS-Specific Recombinant antibody, PBS Only (Detector)

Catalog Number:82682-10-PBS



## Basic Information

Catalog Number:

82682-10-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC006499

GeneID (NCBI):

3265

UNIPROT ID:

P01112

Full Name:

v-Ha-ras Harvey rat sarcoma viral  
oncogene homolog

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

242062A4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01568-1, Human HRAS Recombinant Matched Antibody Pair - PBS only. 82682-11-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00167. 82682-10-PBS was HRP conjugated as the detection antibody. Range: 31.3-2000 pg/mL.



The mean HRAS concentration was determined to be 17.2 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load.