

For Research Use Only

HRAS-Specific Recombinant antibody, PBS Only (Capture)

Catalog Number:82682-11-PBS



Basic Information

Catalog Number: 82682-11-PBS	GenBank Accession Number: BC006499	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3265	CloneNo.: 242062B5
Source: Rabbit	UNIPROT ID: P01112	
Isotype: IgG	Full Name: v-Ha-ras Harvey rat sarcoma viral oncogene homolog	
	Calculated MW: 21 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01568-1, Human HRAS Recombinant Matched Antibody Pair - PBS only. 82682-11-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00167. 82682-10-PBS was HRP conjugated as the detection antibody. Range: 31.3-2000 pg/mL.



The mean HRAS concentration was determined to be 17.2 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load.