

For Research Use Only

MMP9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82685-19-PBS

Basic Information

Catalog Number: 82685-19-PBS	GenBank Accession Number: BC006093	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4318	CloneNo.: 241224E6
Source: Rabbit	UNIPROT ID: P14780	
Isotype: IgG	Full Name: matrix metalloproteinase 9 (gelatinase B, 92kDa gelatinase, 92kDa type IV collagenase)	
	Calculated MW: 707 aa, 78 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01003-1, MMP-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82685-18-PBS. Detection antibody: 82685-19-PBS. Standard: Eg0442. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01003-3, MMP-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82685-17-PBS. Detection antibody: 82685-19-PBS. Standard: Eg0442. Range: 0.625-80 ng/mL