

For Research Use Only

EIF6 Recombinant antibody, PBS Only



Catalog Number: 82693-1-PBS

Basic Information

Catalog Number: 82693-1-PBS	GenBank Accession Number: BC001119	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3692	CloneNo.: 2G23
Source: Rabbit	UNIPROT ID: P56537	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 6	
Immunogen Catalog Number: AG0324	Calculated MW: 27 kDa	
	Observed MW: 27 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

p27(BBP/eIF6) is an evolutionarily conserved protein that was originally identified as p27(BBP), it functions as an interactor of the cytoplasmic domain of integrin 4 and as the putative translation initiation factor eIF6. p27BBP is found in two pools: one nuclear pool enriched in the perinucleolar region, and one cytoplasmic pool. p27BBP binds to the fibronectin type III domains of integrin 4 subunit (ITGB4), an important functional component of hemidesmosomes, and help link ITGB4 to the intermediate filament cytoskeleton. In vitro and in vivo studies demonstrated that p27BBP is essential for cell viability and has a primary function in the biogenesis of the 60S ribosomal subunit. p27BBP protein is increased in rapidly cycling cells and decreased in villous cells committed to apoptotic cell death. In dysplastic colorectal adenomas and carcinomas, p27BBP displayed a large increase of its nucleolar component and was associated with the nuclear matrix. In particular, p27BBP increased progressively from adenomas to carcinomas and was related to the tumor stage.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82693-1-RR (EIF6 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82693-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82693-1-RR against Human EIF6 were performed. The affinity constant is 6.29 nM.