

For Research Use Only

IL1B Recombinant antibody, PBS Only (Capture)



Catalog Number: 82696-3-PBS

Basic Information

Catalog Number: 82696-3-PBS	GenBank Accession Number: BC008678	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3553	CloneNo.: 304
Source: Rabbit	ENSEMBL Gene ID: ENSG00000125538	
Isotype: IgG	UNIPROT ID: P01584	
	Full Name: interleukin 1, beta	
	Calculated MW: 269 aa, 31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80001-3, IL-1 beta Hybrid Matched Antibody Pair - BSA and Azide free. Capture antibody: 82696-3-RR. Detection antibody: 68767-2-PBS. Standard:HZ-1164. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP80001-5, IL-1 beta Hybrid Matched Antibody Pair - BSA and Azide free. Capture antibody: 82696-3-RR. Detection antibody: 65276-1-PBS. Standard:HZ-1164. Range: 3.125-100 ng/mL