

For Research Use Only

IL1B Recombinant antibody, PBS Only (Capture)



Catalog Number: 82696-7-PBS

Basic Information

Catalog Number:

82696-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC008678

GeneID (NCBI):

3553

ENSEMBL Gene ID:

ENSG00000125538

UNIPROT ID:

P01584

Full Name:

interleukin 1, beta

Calculated MW:

269 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

3C16

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00331-1, IL-1 beta Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 82696-7-PBS. Detection antibody: 82696-5-PBS. Standard:HZ-1164. Range: 0.781-100 ng/mL