

For Research Use Only

IL1B Recombinant antibody

Catalog Number: 82696-8-RR



Basic Information

Catalog Number: 82696-8-RR	GenBank Accession Number: BC008678	Purification Method: Protein A purification
Size: N/A	GeneID (NCBI): 3553	CloneNo.: 2G9
Source: Rabbit	ENSEMBL Gene ID: ENSG00000125538	
Isotype: IgG	UNIPROT ID: P01584	
	Full Name: interleukin 1, beta	
	Calculated MW: 269 aa, 31 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data