

For Research Use Only

# Cyclin E Recombinant antibody, PBS Only (Capture)



Catalog Number: 82700-5-PBS

## Basic Information

Catalog Number: 82700-5-PBS	GenBank Accession Number: BC035498	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 898	CloneNo.: 240136C5
Source: Rabbit	UNIPROT ID: P24864	
Isotype: IgG	Full Name: cyclin E1	
Immunogen Catalog Number: AG2110	Calculated MW: 410 aa, 47 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00361-3, Cyclin E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82700-5-PBS. Detection antibody: 82700-8-PBS. Standard: Ag2110. Range: 0.625-40 ng/mL.