

For Research Use Only

# Cleaved Caspase 3; P17

## Recombinant antibody, PBS Only (Capture)

Catalog Number: **82707-8-PBS**

### Basic Information

Catalog Number: 82707-8-PBS	GenBank Accession Number: NM_032991	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 836	CloneNo.: 230527B9
Source: Rabbit	UNIPROT ID: P42574	Affinity: $K_D = 1.47 \times 10^{-9}$
Isotype: IgG	Full Name: caspase 3, apoptosis-related cysteine peptidase	$K_{off} = 9.91 \times 10^{-4}$ $K_{on} = 6.73 \times 10^5$
	Calculated MW: 32 kDa	

### Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

### Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00104-3, Cleaved Caspase 3/P17 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 82707-8-PBS. Detection antibody: 82707-11-PBS. Standard: SY01603. Range: 0.156-2.5 ng/mL