

For Research Use Only

# IL22 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82766-13-PBS

## Basic Information

Catalog Number:

82766-13-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC070261

GeneID (NCBI):

50616

UNIPROT ID:

Q9GZX6

Full Name:

interleukin 22

Calculated MW:

179 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

240699C11

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00657-1, IL-22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82766-12-PBS. Detection antibody: 82766-13-PBS. Standard: Eg0228. Range: 0.313-40 ng/mL