

For Research Use Only

# Mouse Ifng Recombinant antibody, PBS Only (Detector)



Catalog Number: 82768-4-PBS

## Basic Information

Catalog Number:

82768-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC119063

GeneID (NCBI):

15978

UNIPROT ID:

P01580

Full Name:

interferon gamma

Calculated MW:

17 kDa

Purification Method:

Protein A purification

CloneNo.:

240262B6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00498-1, MOUSE Ifng Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82768-6-PBS. Detection antibody: 82768-4-PBS. Standard: Eg0334. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00498-2, MOUSE Ifng Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82768-5-PBS. Detection antibody: 82768-4-PBS. Standard: Eg0334. Range: 0.78-100 ng/mL