

For Research Use Only

# Beta Arrestin 2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82780-7-PBS



## Basic Information

Catalog Number: 82780-7-PBS	GenBank Accession Number: BC007427	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 409	CloneNo.: 241470G2
Source: Rabbit	UNIPROT ID: P32121	
Isotype: IgG	Full Name: arrestin, beta 2	
Immunogen Catalog Number: AG0227	Calculated MW: 46 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01125-4, Human Beta Arrestin 2 Recombinant Matched Antibody Pair - PBS only. 82780-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0227. 82780-4-PBS was HRP conjugated as the detection antibody. Range: 0.313-10 ng/mL



The mean Beta Arrestin 2 concentration was determined to be 9.58 ng/mL in HT-29 cell extract based on a 1.30 mg/mL extract load and 5.35 ng/mL in mouse brain cell extract based on a 2.00 mg/mL extract load.