

For Research Use Only

OX40L/TNFSF4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82794-2-PBS



Basic Information

Catalog Number:

82794-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003326.5

GeneID (NCBI):

7292

UNIPROT ID:

P23510

Full Name:

tumor necrosis factor (ligand)

superfamily, member 4

Purification Method:

Protein A purification

CloneNo.:

240199E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00327-1, OX40L/TNFSF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82794-4-PBS. Detection antibody: 82794-2-PBS. Standard: Eg0582. Range: 0.313-40ng/mL.



Biolayer interferometry (BLI) kinetic assay of 82794-2-PBS against Human OX40L/TNFSF4 was performed. The affinity constant is 0.217 nM.