

For Research Use Only

# HNRNPC Recombinant antibody, PBS Only (Capture)



Catalog Number: 82798-5-PBS

## Basic Information

Catalog Number:

82798-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2356

GenBank Accession Number:

BC003394

GeneID (NCBI):

3183

UNIPROT ID:

P07910

Full Name:

heterogeneous nuclear  
ribonucleoprotein C (C1/C2)

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

230373F6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00088-1, HNRNPC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82798-5-PBS. Detection antibody: 82798-4-PBS. Standard: Ag2356. Range: 1.25-20 ng/mL.