

For Research Use Only

Mouse ICAM1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82827-3-PBS

Basic Information

Catalog Number:

82827-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010493

GeneID (NCBI):

15894

UNIPROT ID:

P13597

Full Name:

intercellular adhesion molecule 1

Calculated MW:

59kd

Purification Method:

Protein A purification

CloneNo.:

240249H3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00543-1, MOUSE ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82827-6-PBS. Detection antibody: 82827-3-PBS. Standard: Eg0482. Range: 0.313-40 ng/mL.