

For Research Use Only

BACH1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82842-1-PBS

Basic Information

Catalog Number:

82842-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5128

GenBank Accession Number:

BC063307

GeneID (NCBI):

571

ENSEMBL Gene ID:

ENSG00000156273

UNIPROT ID:

O14867

Full Name:

BTB and CNC homology 1, basic
leucine zipper transcription factor 1

Calculated MW:

82 kDa

Observed MW:

100-110 kDa

Purification Method:

Protein A purification

CloneNo.:

4P21

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80007-2, BACH1 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66762-2-PBS. Detection antibody: 82842-1-PBS. Standard: Ag6011. Range: 0.313-40 ng/mL