

For Research Use Only

Mouse Snf1lk Recombinant antibody, PBS Only



Catalog Number: 82851-3-PBS

Basic Information

Catalog Number: 82851-3-PBS	GenBank Accession Number: NM_010831	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 17691	CloneNo.: 3D2
Source: Rabbit	UNIPROT ID: Q60670	
Isotype: IgG	Full Name: salt inducible kinase 1	
Immunogen Catalog Number: AG33435	Calculated MW: 85 kDa	
	Observed MW: 85 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82851-3-RR (Snf1lk antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82851-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82851-3-RR against Mouse Snf1lk were performed. The affinity constant is 0.686 nM.