

For Research Use Only

TNFRSF17 Recombinant antibody, PBS Only



Catalog Number: 82855-1-PBS

Basic Information

Catalog Number: 82855-1-PBS	GenBank Accession Number: BC058291	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 608	CloneNo.: 5H7
Source: Rabbit	UNIPROT ID: Q02223	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 17	
	Calculated MW: 184 aa, 20 kDa	
	Observed MW: 20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82855-1-RR (TNFRSF17 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82855-1-PBS in a different storage buffer formulation.