

For Research Use Only

# BCMA/TNFRSF17 Recombinant antibody, PBS Only

Catalog Number:82855-1-PBS



## Basic Information

Catalog Number:

82855-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC058291

GeneID (NCBI):

608

UNIPROT ID:

Q02223

Full Name:

tumor necrosis factor receptor  
superfamily, member 17

Calculated MW:

184 aa, 20 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

5H7

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82855-1-RR (TNFRSF17 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82855-1-PBS in a different storage buffer formulation.