

For Research Use Only

TLN1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82856-14-PBS

Basic Information

Catalog Number: 82856-14-PBS	GenBank Accession Number: BC042923	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7094	CloneNo.: 240294F10
Source: Rabbit	UNIPROT ID: Q9Y490	
Isotype: IgG	Full Name: talin 1	
Immunogen Catalog Number: AG5367	Calculated MW: 2541 aa, 270 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00405-1, Talin-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82856-13-PBS. Detection antibody: 82856-14-PBS. Standard: Ag5367. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00405-2, Talin-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82856-14-PBS. Detection antibody: 82856-17-PBS. Standard: Ag5367. Range: 1.56-100 ng/mL