

For Research Use Only

# TLN1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82856-17-PBS

## Basic Information

Catalog Number: 82856-17-PBS	GenBank Accession Number: BC042923	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7094	CloneNo.: 240294F5
Source: Rabbit	UNIPROT ID: Q9Y490	
Isotype: IgG	Full Name: talin 1	
Immunogen Catalog Number: AG5367	Calculated MW: 2541 aa, 270 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00405-2, Talin-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82856-14-PBS. Detection antibody: 82856-17-PBS. Standard: Ag5367. Range: 1.56-100 ng/mL