

For Research Use Only

Mouse Gm-csf Recombinant antibody, PBS Only (Detector)



Catalog Number:82861-7-PBS

Basic Information

Catalog Number:

82861-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM-009969

GeneID (NCBI):

12981

UNIPROT ID:

P01587

Full Name:

colony stimulating factor 2
(granulocyte-macrophage)

Purification Method:

Protein A purification

CloneNo.:

230286A6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00306-1, Mouse GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82861-9-PBS. Detection antibody: 82861-7-PBS. Standard: Eg0512. Range: 0.625-10 ng/mL.