

For Research Use Only

Mouse Gm-csf Recombinant antibody, PBS Only (Capture)



Catalog Number: 82861-9-PBS

Basic Information

Catalog Number: 82861-9-PBS	GenBank Accession Number: NM-009969	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 12981	CloneNo.: 230286E11
Source: Rabbit	UNIPROT ID: P01587	
Isotype: IgG	Full Name: colony stimulating factor 2 (granulocyte-macrophage)	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00306-1, Mouse GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82861-9-PBS. Detection antibody: 82861-7-PBS. Standard: Eg0512. Range: 0.625-10 ng/mL.



Cytometric bead array standard curve of MP00306-2, Mouse GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82861-9-PBS. Detection antibody: 82861-8-PBS. Standard: Eg0512. Range: 0.625-10 ng/mL.