

For Research Use Only

# Mouse Il12b Recombinant antibody, PBS Only (Detector)



Catalog Number:82865-6-PBS

## Basic Information

Catalog Number:

82865-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001303244

GeneID (NCBI):

16160

UNIPROT ID:

P43432

Full Name:

interleukin 12b

Calculated MW:

38kd

Purification Method:

Protein A purification

CloneNo.:

240210E4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00413-1, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-8-PBS. Detection antibody: 82865-6-PBS. Standard: Eg0299. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00413-3, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-9-PBS. Detection antibody: 82865-6-PBS. Standard: Eg0299. Range: 5-80 ng/mL.