

For Research Use Only

UGT8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82874-4-PBS

Basic Information

Catalog Number: 82874-4-PBS	GenBank Accession Number: BC075069	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7368	CloneNo.: 230107G8
Source: Rabbit	UNIPROT ID: Q16880	
Isotype: IgG	Full Name: UDP glycosyltransferase 8	
Immunogen Catalog Number: AG13137	Calculated MW: 541 aa, 61 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00003-1, UGT8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82874-4-PBS. Detection antibody: 82874-3-PBS. Standard: Ag13137. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00003-2, UGT8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82874-4-PBS. Detection antibody: 82874-5-PBS. Standard: Ag13137. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00003-3, UGT8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82874-4-PBS. Detection antibody: 82874-2-PBS. Standard: Ag13137. Range: 0.625-20 ng/mL